

PROTODERM™ XEROFORM PETROLATUM PIN SITE DRESSINGS

The management of skeletal pin sites (for patients with external fixation devices) represents a major dilemma for orthopaedic nurses all over the world. Infection rates are high and wound management costs are significant.

Integrity Medical Devices, Inc. is proud to introduce ProtoDerm® Pin Site Dressings, the only Xeroform Petrolatum dressing available specifically designed to help manage external fixation device pin site wounds. PROTODERM is used as a primary layer to cover and protect wounds from environmental contamination with mild deodorizing.

The unique die-cut individual dressings, allows for easy access and precise pin/wire placement. PROTODERM is the only product packaged 10 ready-to-use dressings in one, sterile pouch resulting in significantly lower cost per dressing and ease of use. The dressings are easily pulled apart from the single packaged strip of 10. Each sterile pouch contains 10, 1" x 1" square ProtoDerm Pin Site Dressings. Each dressing is slit and slotted to fit easily around external fixation pins and wires.



Now fixator pin care is as easy as *One, Two, Three.*


One – Peel Open Sterile Package

Each package contains 10 individual 1" x 1" pin dressings.



Two – Pull apart a 1" x 1" ProtoDerm Dressing

Each dressing pulls apart easily from the main pad of 10 dressings.



Three – Slide Dressing around Pin

The custom slit and pin slot slide easily onto the pin fit

